


Non Contact Infrared Body
THERMOMETER

QM SMART CLOUD CO., LIMITED

Measurement Method



Measurement Method

Voltage	DC3V,AAA battery (1.5V * 2)		
Measuring range	Body temperature mode 34-42.9°C		
Measurement accuracy	34.0-34.9°C:±0.3°C 35.0-42.0°C:±0.2°C 42.1-42.9°C:±0.3°C		
Measuring distance	1-3 cm	Measure time	≤2 seconds
Shut down	≤18 seconds	Memorize data	32 groups
Storage and Transport Conditions	-20 °C to 55 °C at a relative humidity of≤ 93%		
Operating Conditions	-16°C-35°C at a relative humidity of≤85%		

